

# Tourniquet

## Range of application


For the vascular hemorrhage on wounded limbs, the battle line emergency firstaid, also used in pre-hospital emergency treatment.

## Specification

Unit:mm

Specification	basic sizes			
	width	Tolerances	length	Tolerances
38X950mm	38	$\pm 2$	950	$\pm 10$
38X750mm	38	$\pm 2$	750	$\pm 10$
38X650mm	38	$\pm 2$	650	$\pm 10$

## Operating procedures



A — Rub Buckle

B — Rod

D — Adhesive Tape


C — Fixed Buckle

E — Fixed Slot

### Arm bleeding one-handed operation

- ① Take it out of the box and open the adhesive tape, sheathe it to the wound closed to the heart
- ② Tighten the adhesive tape, stick tight reversely
- ③ Stuck tight along the upper arm, don't cover the screw rods
- ④ Rotating the rod, until stopped bleeding
- ⑤ Fix the rod in slot, then the bleeding have been controlled
- ⑥ As to the upper limb or calf, continue glue tightly and placed on the rod
- ⑦ Lock the screw rods and adhesive tape fixed clasp

1



2



3



4



5



6




7





## Double-hand operation legs bleeding

- ① When arounding the legs ,let the free end through the inside buttonhole
- ② Back cross the tape through outer buttonhole ,fx the band with friction
- ③ Tighten the adhesive tape ,reversely tighted
- ④ Rotating the rod ,until the bleeding stopped
- ⑤ Stuck the rod with the fixed slot, then the bleeding have been controlled
- ⑥ Lock the screw rods with the fixed clasp



## Usage of friction buckle

When manually operated a tourniquet or when the tourniquet is serious polluted, use special design of friction belt buckle to fix the band.

- a. Let the adhesive tape free end through buckles inside
- b. Back through the outer buttonhole, so the adhesive tape can be firmly positioning
- c. Tighten the adhesive tape, reversely tighted

## Cautions

- a. Disposables, single use.
- b. This tourniquet be used by medical first-aid trained personnel.
- c. Read the instructions before using.
- d. Medical waste disposal requirements: according to local regulations.
- e. According to medical wastes treatment, not discarded.

## Storage

- a. Under the storage conditions, the period of validity is 5 years.
- b. Generally should be stored in a temperature-5- 30 degrees Celsius, relative humidity no more than 80%, non- corrosive gases and ventilated interior.